


#### Introduction to Single Particle Reconstruction

#### Steve Ludtke


Charles C. Bell Professor
Biochemistry and Molecular Biology
Director, CryoEM/CryoET Core
Co-director CIBR Center
Baylor College of Medicine







#### 190/~30,000 Particle Images

Too, co,coo i article illiages																	
P		4		+		1 2 5						ů.		400			
											7				(A)		
				2					19.E	200	×						
*		1.															
					V.									300			c.
																	15
				3.					er.	7,4				*			
77.						4-1	***			i.e						9.	
					1												46
								法					7.35				107



Bai, X. C., Fernandez, I. S., McMullan, G. & Scheres, S. H. Ribosome structures to near-atomic resolution from thirty thousand cryo-EM particles. Elife 2, e00461 (2013). PMC3576727.






S. Cereviseae 80S Ribosome (EMD-2275)


Data taken from PDBe 3DEM test data
no movie alignment performed

Dataset 10002 (Bai XC, Fernandez IS,

McMullen G, Scheres SH)




#### TRPV1



Liao, M., Cao, E., Julius, D., and Cheng, Y. (2013). Structure of the TRPV1 ion channel determined by electron cryo-microscopy. Nature. 504:107-112.


# 30 kDa HIV-1 RNA Dimerization Signal





~9 Å

# 30 kDa HIV-1 RNA Dimerization Signal










~9 Å EMAN2

				7				
								(V)
			120					
Fatty	Acid Sy	nthase						


*	*		<b>33</b>	3	8	<b>CR</b>
33	1052		977	0	•	2
	722	748	910	045	6	220
778	4	<b>670</b>		200		
<b>S</b>	633	1199	1390	9	1201	(3)


## Unsupervised Classification




GroEL single ring mutant + GroES

Fatty Acid Synthase ~30 Å motion










IP3 Receptor ~18 Å motion




TriC (open) ~12 Å motion












#### β Galactosidase








β Galactosidase

(this is real ~2.2 Å B-gal data)

## Unsupervised Classification




GroEL single ring mutant + GroES


### Invariants?

- An image derived from another image which doesn't change under some operation
- A translational invariant is an image that doesn't change when its source is moved around in the box
- A rotational invariant is an image that doesn't change when its source is rotated


If we had a good one we could skip costly alignments!

## Autocorrelation







# Autocorrelation



β Galactosidase

# Autocorrelation









Any One of These


And we still haven't done rotation....

# Problems with Autocorrelation

- We lose phases, twice. Significant information loss
- Handedness lost
- Noise additivity
  - P(A+B+C+D) != P(A)+P(B)+P(C)+P(D)







(this is real ~2.2 Å B-gal data)

· Power spectrum:

$$P(\bar{s}) = |F(\bar{s})|^2$$

Bispectrum

$$B(\bar{s}_1, \bar{s}_2) = F(\bar{s}_1)F(\bar{s}_2)F^*(\bar{s}_1 + \bar{s}_2)$$
  
**4-D**!?!

3 point correlation in real space translationally invariant

Power spectrum:


$$P(\bar{s}) = |F(\bar{s})|^2$$

Bispectrum

$$B(\bar{s}_1, \bar{s}_2) = F(\bar{s}_1)F(\bar{s}_2)F^*(\bar{s}_1 + \bar{s}_2)$$
**4-D** !?!

3 point correlation in real space


... So we compute a carefully selected 3-D subspace from the bispectrum integrated rotationally.



Significant changes with 3-D orientation



β Galactosidase



Good translational invariance




β Galactosidase



Good rotational invariance





Good noise behavior




β Galactosidase


# Bispectral Invariants

- Preserve more information
- Good noise resistance
- Noise is still noise and can cancel out with averaging
- Relatively inexpensive (generate ~50 each sec per CPU)


## (old) 2D Refinement



# (new) 2D Refinement




TcdA1




5513 particles 240x240 box size 64 Classes

2 h 27 m 1 Laptop (4 core)




5513 particles 240x240 box size 64 Classes

7 min 1 Laptop (4 core)




8,290 particles 384x384 box size 128 Classes 15 sec to classify 4 min total 1 workstation



10,533 particles 208x208 box size 128 Classes

2h 39m 1 laptop



10,533 particles 208x208 box size 128 Classes

12 min 1 laptop

## Single Particle Analysis

Prepare Grids Record Images Movie Alignment **Extract Particles** Determine Orientations 3-D Reconstruction Interpret Map

## Single Particle Analysis

**Extract Particles** 



Determine Orientations





3-D

Reconstruction




Interpret Map





7.5° Angular Step

#### **Determine Particle Orientations**



Every particle vs every possible projection?

#### Determine Particle Orientations

(Ribosome at ~4 Å)

- \* 30,000 particles
- \* 3000 projections (~2.6 degree sampling)
- \* 180 rotations
- \* 10 x 10 translations (if centering is decent)
- \* 65,536 pixels (256x256 image)
- \* 20 FLOPS/pixel
- = 2.1 x 1018 FLOPS (2 exaFLOPS)
- @ 100 gigaFLOPS/s = 6,000 hours (~ 240 days on one typical workstation)
- Clearly we need to be smarter...
- (Actual EMAN2 time <3 hours on ... even more cleverness possible though)</li>

#### **Determine Particle Orientations**

- \* 30,000 particles
- \* 3000 projections (~2.6 degree sampling)
- \* 180 rotations
- \* 10 x 10 translations (if centering is decent)
- \* 65,536 pixels (256x256 image)
- \* 20 FLOPS/pixel
- = 2.1 x 1018 FLOPS (2 exaFLOPS)
- @ 100 gigaFLOPS/s = 6,000 hours (~ 240 days on one typical workstation)
- Clearly we need to be smarter...
- (Actual EMAN2 time <3 hours ... even more cleverness possible though)</li>

## Single Particle Analysis

**Extract Particles** 



Determine Orientations



3-D

Reconstruction





Interpret Map

# Reconstruction Algorithms


- Back Projection
- Filtered Back Projection
- Direct Fourier Inversion
- SIRT
- SART


• ...








Crystals have spots
Particles are ~continuous and have phases!











# One problem...

- To determine the particle orientations, we needed a 3-D reference
- catch 22
- All current CryoEM Single Particle solutions are iterative!